

## Resolving Buffer

产品货号: R22039

产品规格: 500ml

### 产品简介:

Resolving Buffer Packs are pouches of dry blend powder each sufficient to make 500 mL of resolving gel buffer (1.5 M Tris buffer, pH 8.8) for handcasting polyacrylamide gels. They are easy to use—simply empty the contents of one pouch into a beaker, add ultrapure water, and stir to dissolve.

Features of the SureCast Resolving Buffer Packs include:

- Convenient—dissolve contents of one packet in water and the buffer is ready to use
- Saves time and space—no weighing, no calculations, no pH adjustment, no need to stock individual components
- Long shelf life—stocking and storage as dry powder eliminates concerns about long-term stability of stock solutions

### Applications

Resolving gel buffer for casting polyacrylamide gels



扫一扫 加微信

上海尚宝生物科技有限公司

Shanghai Saint-Bio Biotechnology Co., Ltd

地址: 上海市徐汇区龙华路2518弄14号

免费电话: 400-611-0007 13671551480 13643719779

Q Q: 807961520 731791866

邮箱: shsunbao@126.com

<http://www.saint-bio.com>