

## 20× Borate Buffer

产品货号: R21918

产品规格: 100ml/500ml

### 产品简介:

20× Borate Buffer is a space-saving stock solution that is ideal for preparing sodium borate buffer solutions for use in protein modification procedures requiring amine-free buffer at alkaline pH.

Concentrated Buffers are ready to use without having to weigh and dissolve dry ingredients or to adjust the pH with concentrated acid or base. Simply dilute the stock solution with pure water and proceed with your experiment. The 20× Borate Buffer makes 50mM borate, pH 8.5, when diluted to 1× with water.

### 产品组成:

产品名称	规格	保存条件
20× Borate Buffer	100ml/500ml	室温

### 产品特点:

- Borate buffer—diluted 20-fold in water, the solution yields 50mM borate, pH 8.5.
- Easy to use—no packets to open and no powder to dissolve.
- Increased accuracy—eliminates the possibility of powder remaining in a packet.
- Saves time—20× concentration eliminates time spent waiting for powder to dissolve.
- Saves space—storage as concentrated stock minimizes bench space needed for solutions.

### 产品应用:

- Protein modification procedures requiring amine-free buffer at alkaline pH.



扫一扫 加微信

上海尚宝生物科技有限公司

Shanghai Saint-Bio Biotechnology Co.,Ltd

地址:上海市徐汇区龙华路2518弄14号

电话:400-611-0007 13671551480

Q Q: 807961520

邮箱: saintbio@126.com

<http://www.saint-bio.com>